




Bürettenspitze / Schlauchlänge 0.4 m

6.1537.010

Lieferumfang 6.1537.010

Qt.	Order no.	Beschreibung
1 PCS	6.1446.030	Ball stopper / B-14/(15) For the movable holding of buret tips in titration vessels



1 PCS

6.1543.060

Tip / M6 thread

Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.



1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection

