



Probenrack, 48 x 11 mm, (Al)

6.2041.200

Probenrack zu 889 IC Sample Center passend für 48 Probengefäße. Die folgenden Probengefäße können verwendet werden: 6.2743.100, 6.2743.107, 6.2743.110, 6.2743.117, 6.7303.100.

Material: Aluminium (Al)

Optionales Zubehör

Order no.	Beschreibung
6.2743.100	PP Vials 0.7 mL to IC Sample Center, 1000 pieces 1000 sample vials.



6.2743.107 PP Vials 0.7 mL to IC Sample Center, 100 pieces
100 sample vials.




6.2743.110 PP Vials 0.3 mL to IC Sample Center, 1000 pieces
1000 sample vials.



6.2743.117 PP Vials 0.3 mL to IC Sample Center, 100 pieces
100 sample vials.



6.2743.120 Septum to 6.2743.1XX Si/PTFE
1000 pieces



6.2743.127

Septum to 6.2743.1XX Si/PTFE

100 pieces



6.7303.100

Glass vial, 1.5 mL, 100 pieces

For use with the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ) and 889 Sample Center.



6.7303.110

Septum for 67303100, 100 pieces

For use with the 67303100 Glass vial on the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ), and 889 Sample Center.

