



Pufferlösung pH 7 (500 mL)

6.2307.110

gebrauchsfertige Pufferlösung pH = 7 (25 °C), eingefärbt, mit Verschlussiegel

Optionales Zubehör

Order no.	Beschreibung
6.2307.100	Buffer solution pH 4 (500 mL) Ready-to-use buffer solution pH = 4 (25 °C), colored, with closure seal



6.2307.120 Buffer solution pH 9 (500 mL)
Ready-to-use buffer solution pH = 9 (25 °C), colorless, with closure seal



6.2307.200 Buffer solution pH 4.00 (30 x 30 mL)
Buffer solution pH 4.00 (25 °C) in single use sachets, colorless



6.2307.210 Buffer solution pH 7.00 (30 x 30 mL)
Buffer solution pH 7.00 (25 °C) in single use sachets, colorless



6.2307.220 Buffer solution pH 9.00 (30 x 30 mL)
Buffer solution pH 9.00 (25 °C) in single use sachets, colorless



6.2307.230

Buffer solutions pH 4, 7 and 9

Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL



6.2325.000

pHit kit

Maintenance kit for pH electrodes

The kit contains:

- 50 mL cleaning solution
- 50 mL 3M KCl solution
- 50 mL storage solution
- 2 Storage vessels
- Instructions for use

