



Austauschset O₂-Lumitrode

6.5623.000

Austauschset für die O₂-Lumitrode. Das Set beinhaltet eine zertifizierte O₂-Kappe als Ersatz, einen Kalibrierköcher und 5 x 30 mL Sauerstoffstandard, 0%.

Lieferumfang 6.5623.000

Qt.	Order no.	Beschreibung
1 PCS	6.1260.000	Calibration vessel O2
		Calibration vessel for calibrating the O ₂ -Lumitrode at 100% air saturation.



Oxygen standard for 0% calibration of the O₂-Lumitrode. Contains 5 sachets of 30 mL each.



Optionales Zubehör

Order no.	Beschreibung
6.1116.000	O2 Lumitrode The optical sensor for measuring dissolved oxygen (DO) can be used with a 913 pH/DO meter or with a 914 pH/DO conductometer. The measuring principle of the sensor is based on luminescence quenching. The space-saving and maintenance-free sensor is suitable for DO measurement, for example in: • Water quality control • Wastewater industry • Beverage production • Fish farming



This sensor is supplied with a calibration vessel.

Where necessary, it is easy to replace the measurement cap (O₂ cap) that contains the oxygen-intensive luminophore.
