


BWID Standard-Raman-ID-Software

BWT-840000126


Software zur Datenaufnahme, Materialidentifikation und Erstellung von Reports für die tragbaren Raman-Spektrometer von B&W Tek, die eine schnelle Identifikation oder Verifizierung von Materialien auf Grundlage von benutzererstellten Bibliotheken, urheberrechtlich geschützten Bibliotheken von B&W Tek sowie Bibliotheken von Drittanbietern ermöglicht. Inklusive Leistungsüberprüfung des Systems. Umfasst eine Benutzerlizenz. Kostenlose Upgrades sind für ein Jahr ab Kaufdatum inbegriffen.

Optionales Zubehör

Order no.	Beschreibung
BWT-840001123	Bluetooth 1D barcode scanner
	Bluetooth 1D barcode scanner for use with BWID software operating a portable Raman instrument



BWID library of illicit materials includes drugs, TICs, TIMs, precursors, and general chemicals.



BWID Metrohm Comprehensive RamanLibrary - More than 27,700 spectral entries of chemicals, minerals, organic chemicals, solvents, polymers, dyes and pigments, TICs, TIMs, inorganics and organometallics, personal care products, and illicit materials.

PTL-MCRL

