



TX-900MH impact printer

2.143.0200

Compact impact printer with USB and RS232 interface for

- 900 Touch Control (+)
- 915 KF Ti-Touch (+)
- 916 Ti-Touch (+)
- 917 Coulometer (+)
- 91X Meter (+)
- Eco Dosimat/Titrator (+)
- 877/848 Titrino plus (additional 6.2151.100 required)
- 865/876 Dosimat plus (additional 6.2151.100 required)
- 862 Compact Titrosampler (additional 6.2151.100 required)
- 870 KF Titrino plus (additional 6.2151.100 required)
- 899 Coulometer (additional 6.2151.100 required)

(+) Cable included in scope of delivery

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.143.0200

Qt.	Order no.	Description
1 PCS	1.143.0200	TX-900MH impact printer
		Compact impact printer with USB and RS232 interface for
		900 Touch Control
		915 KF Ti-Touch
		916 Ti-Touch
		917 Coulometer
		877/848 Titrino plus
		865/876 Dosimat plus
		91X meter (cable 6.2151.140)
		Eco Dosimat
		Eco Titrator
		Eco Coulometer
		862 Compact Titr sampler
		870 KF Titrino plus
		899 Coulometer
1 PCS	6.2134.040	RS-232 DB9 - IBM-PC(DB9) connecting cable
		For the Metrohm RS-232 interface - PC connection. For Titrinos, Titrino plus and 756, 831 KF Coulometer, RS-232/USB Box
1 PCS	6.2151.020	Cable USB A - USB B 1.8 m
		USB connecting cable



Y cable for connection of a USB printer to the pH/Conductometers 912 / 913 / 914. This cable allows a printer and the power supply unit to be connected to the measuring instrument at the same time.



Optional accessories

Order no.	Description	
6.2151.100	Adapter USB MINI (OTG) - USB A For connecting USB instruments.	
6.9001.070	Paper roll for TX-900MH Paper roll for TX-900MH impact printer	
6.9001.080	Continuous adhesive labels for TX-900MH Self-adhesive paper roll for TX-900MH impact printer	
6.9002.070	Color ribbon for TX-900MH impact printer	