



700 Dosino

2.700.0010

Drive for dosing unit, with permanently installed cable and 9-pin plug (DB9).

This product is no longer available!

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.700.0010

Qt.	Order no.	Description
1 PCS	1.700.0010	700 Dosino
		Drive for Dosing Unit, with fixed cable and 9-pin plug (DB9).



1 PCS

6.1546.030

Piston tongs

For pistons for 807 Dosing Units or cylinder units OMNIS.



1 PCS

6.2134.200

Adapter DB9 St. /Mini Din8 Bu



1 PCS

6.2148.040

Adapter for Dosino 700 (series 23 or higher)



Optional accessories

Order no.	Description	
6.2134.020	Adapter DB9 female/mini-DIN8 male DB9 socket and mini-DIN 8 plug. Connection 700.0010 Dosino to 736, 751, 758, 799 Titrinos, 756, 831 KF Coulometer, 711 Liquino.	