



# 700 Dosino

2.700.0020

The 700 Dosino is a drive for dosing units. With a permanently installed cable and Mini DIN plug.

This product is no longer available!

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 2.700.0020

Qt.	Order no.	Description
1 PCS	1.700.0020	700 Dosino
		The 700 Dosino is a Dosing Unit drive. With fixed cable and mini-DIN plug.



1 PCS

6.1546.030

Piston tongs

For pistons for 807 Dosing Units or cylinder units OMNIS.



---


1 PCS

6.2148.040

Adapter for Dosino 700 (series 23 or higher)



# Optional accessories

Order no.	Description	
6.2134.010	Adapter DB9 male / Mini-DIN8 female  DB9 plug and mini-DIN 8 socket. E.g. for the connection of 700.0020 Dosinos with mini-DIN 8 plugs to 726, 796 Titroprocessors.	

---