


791 IC VA Detector

2.791.0010

D

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.791.0010

Qt.	Order no.	Description
1 PCS	1.791.0010	791 IC VA Detector
	D	
		
1 PCS	6.2813.020	Dummy Cell for VA Detector
		Electrical simulation of a voltammetry cell, used for simulation and control purposes. For VA Detectors
		

Optional accessories

Order no.	Description	
6.1204.130	Silver Electrode Tip Metal electrode tip, silver (Ag). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.	
6.1204.140	Gold Electrode Tip Metal electrode tip, gold (Au). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.	