



800 Dosino

2.800.0020

Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 0.65 m).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.800.0020

| Qt. | Order no. | Description |
|-------|------------|----------------------------------|
| 1 PCS | 1.800.0040 | 800 Dosino (cable length 0.65 m) |

Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 0.65 m).



1 PCS

6.1546.030

Piston tongs

For pistons for 807 Dosing Units or cylinder units OMNIS.

