



802 Stirrer for sample changer

2.802.0020

Rod stirrer for sample changer and Sample Processor.

With 104 mm 6.1909.020 stirring propeller and permanently installed cable.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.802.0020

Qt.	Order no.	Description
1 PCS	1.802.0010	802 Rod Stirrer for 727 Ti Stand
		Rod stirrer with separate 6.2137.000 Speed controller and 6.1909.010 Propeller stirrer.



1 PCS

6.1909.020

Stirring propeller / 95 mm

Propeller stirrer, length from lower edge of SGJ: 104 mm. For use with beakers with 722, 802 Rod Stirrer



Optional accessories

Order no.	Description	
6.1909.050	Stirring propeller / 102 mm Propeller for 802 Stirrer when this is used with sample racks for 120 mL beakers (6.1459.300)	
6.1909.060	Stirring propeller / 95 mm Stirring propeller with 3 wings for intensive stirring	