



# 823 Membrane Pump Unit «rinse»

2.823.0030

Membrane pump for direct connection to Sample Processors. The «Rinsing» version contains numerous rinsing accessories.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery 2.823.0030

Qt.	Order no.	Description
1 PCS	1.823.0010	823 Membrane Pump Unit
		Membrane pump for direct connection to Sample Processors.



1 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



---

2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



---

3 PCS

6.1805.060

FEP tubing / M6 / 60 cm

With light and kink protection



---

1 PCS

6.1808.170

Distributor for Sample Processors for rinsing and aspirating

Distributor for installation on Sample Processors for rinsing and aspirating when using externally connected pumps.



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



---

1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



---

3 PCS

6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.

