

824 Easy Sample Changer

2.824.0020



The economic sample changer for routine use. Easy Sample Changer Micro: sample changer for titrations with smaller sample volumes. Up to 24 samples with a maximum volume of 75 mL can be processed in a single run. With the KFT equipment the Easy Sample Changer can be converted to a Karl Fischer sample changer.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 2.824.0020

Qt.	Order no.	Description
1 PCS	1.741.0010	741 Magnetic Stirrer
		Magnetic stirrer for sample changers.
		The magnetic stirrer is located beneath the sample rack.
		
24 PCS	6.1432.210	Sample beaker, clear glass, 75 mL, 1 piece
		Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece
		

1 PCS

6.1453.250

Sample beaker, plastic (PP), 250 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece



1 PCS

6.1458.020

Titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



1 PCS

6.1805.110

FEP tubing / M6 / 80 cm

With light and kink protection



24 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 PCS

6.2034.020

Mounting bracket for Magnetic Stirrer 741

To attach the magnetic stirrer to the «824 Easy Sample Changer



1 PCS

6.2041.340

Sample rack, 24 x 75 mL, (PVC)

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 24 sample beakers. The following sample beakers can be used: 6.1432.210, 6.1459.400.

Plastic: Polyvinyl chloride (PVC)



1 PCS

6.2621.140

Hexagon key 2.5 mm




1 PCS

6.2751.080

Safety shield

Safety shield for single-tower Sample Processors without 786 Swing Head.



1 PCS

6.2752.010

Plug cover

Used with the 717, 730, 760 Sample Changers, IC Sample Processors,
774 Oven Sample Processor



Optional accessories

Order no.	Description
6.1453.220	Sample beaker, plastic (PP), 200 mL, 1 piece Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece

