



# 826 pH mobile with Primatrode

2.826.0020

Portable pH meter with IrDa interface, including carrying strap, batteries and 6.0228.020 Primatrode with NTC (IP67).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 2.826.0020

Qt.	Order no.	Description
1 PCS	1.826.0010	826 pH mobile
		Portable pH meter with IrDa interface, including carrying strap and batteries.



1 PCS

6.2050.000

Carrying strap for pH Meter

Carrying strap for 704, 826 pH meter.



1 PCS

6.2133.000

Batteries LR6, 1.5V, 4 pieces

1.5 V LR6 type. Set of 4 pieces. Used with 704 pH Meter, 773 pH/mV Simulator




1 PCS

6.2621.140

Hexagon key 2.5 mm



## Optional accessories

Order no.	Description	
6.2307.230	Buffer solutions pH 4, 7 and 9  Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL	
6.2325.000	pHit kit  Maintenance kit for pH electrodes  The kit contains: <ul style="list-style-type: none"><li>• 50 mL cleaning solution</li><li>• 50 mL 3M KCl solution</li><li>• 50 mL storage solution</li><li>• 2 Storage vessels</li><li>• Instructions for use</li></ul>	