




pH Lab Irda 240 V AUS

2.827.0017

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.827.0017

Qt.	Order no.	Description	
1 PCS	6.2013.010	Clamping ring For support rods with a diameter of 10 mm.	
1 PCS	6.2016.050	Support rod / 300 mm	
1 PCS	6.2021.020	Electrode holder Electrode holder for 4 electrodes and 2 buret tips	

1 PCS

6.2161.040

Plug-in power supply unit AUS 240V/6V DC

Power supply for 827 pH lab



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

