



# 827 pH lab; 230 V, EU with Primatrode

2.827.0114

Laboratory pH meter with IrDa interface and 6.0228.010 Primatrode

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 2.827.0114

Qt.	Order no.	Description
1 PCS	6.2013.010	Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.050

Support rod / 300 mm



---

1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



---

1 PCS

6.2161.010

Plug-in power supply unit EU 230V /6V DC

Power supply for 827 pH lab



---


1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



## Optional accessories

Order no.	Description	
6.2307.230	Buffer solutions pH 4, 7 and 9  Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL	
6.2325.000	pHit kit  Maintenance kit for pH electrodes  The kit contains: <ul style="list-style-type: none"><li>• 50 mL cleaning solution</li><li>• 50 mL 3M KCl solution</li><li>• 50 mL storage solution</li><li>• 2 Storage vessels</li><li>• Instructions for use</li></ul>	