



MetroSil RP 3 Guard/4.0

6.01070.500

The MetroSil RP 3 Guard/4.0 is used to protect the MetroSil RP 3 - 150/4.0 against contaminations from particles and bacteria.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.02821.010	Cartridge holder for MetroSil RP 3 Guard/4.0 cartridge
	Holder for the MetroSil RP 3 Guard/4.0 cartridge for connection to the analytical separation column.

