



# Beaker adapter for OMNIS 120 mL sample beaker

6.01404.030

Beaker adapter for centering the 120 mL sample beaker in the slide of the Pick&Place modules

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Optional accessories

Order no.	Description
6.01400.200	Sample beakers, plastic (PP), 120 mL, 20 pieces
	Sample beakers, plastic (polypropylene, PP), 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.



6.01400.300 Sample beakers, clear glass, 120 mL, 20 pieces

Sample beakers, clear glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.



---

6.01400.303 Sample beakers, amber glass, 120 mL, 20 pieces

Sample beakers, amber glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces

