

Beaker adapter for OMNIS 75 mL sample beaker

6.01404.040


Beaker adapter for centering the 75 mL sample beaker in the slide of the Pick&Place modules

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description	
6.01402.000	Sample beakers, clear glass, 75 mL, 30 pieces	
	Sample beakers, clear glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces	

6.01402.003 Sample beakers, amber glass, 75 mL, 30 pieces

Sample beakers, amber glass, 75 mL, suitable for 602041040 OMNIS
sample rack, 30 pieces

