



Eco titration vessel lid for 601406XXX titration vessels for Polytron

6.01405.060

Titration vessel lid for volumetric Karl Fischer titration with an Eco KF titrator with Polytron

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01405.060

Qt.	Order no.	Description
1 PCS	6.1454.130	Form seal made of FPM for 61414xxx
		for KF and VA instruments



Optional accessories

Order no.	Description	
6.01408.010	Set of seals for Karl Fischer titration cell Complete set of seals, suitable for 6.05609.000 and 6.05609.010.	
6.02709.010	Stopper Plastic stoppers, suitable for 2.1001.4220 OMNIS Karl Fischer Titrator.	
6.02709.020	Septum stopper Plastic septum stopper, suitable for 6014050X0 titration vessel lids.	
6.02709.030	Screw nipple Plastic screw nipple, suitable for 6014050X0 titration vessel lids.	

6.02709.040

Stopper

Stopper for titration vessel lids 6014050X0.



6.02709.050

Guide sleeve

For fixation of the Polytron in the 601405040 titration vessel upper part.

