



50 – 150 mL titration vessel

6.01406.250

50 – 150 mL titration vessel with screw top for use with 6014050X0 titration vessel lids.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.01405.000	5 x SGJ 14 titration vessel lid OMNIS
	Titration vessel lid with five SGJ 14 openings for 601406XX0 titration vessels with screw cap.



6.01405.010 OMNIS titration vessel lid for volumetric Karl Fischer titration

Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator.



6.01405.040 OMNIS titration vessel lid for 601406XXX titration vessels for Polytron

Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator with Polytron.

