



Cylinder unit OMNIS special 20 mL, without accessories

6.01508.220

Intelligent cylinder unit 20 mL for one OMNIS Titrator, Titration Module or Dosing Module. This cylinder unit is especially recommended for the following solutions:

- Aqueous alkaline solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkaline solutions
- Permanganate solutions
- EDTA solutions

Without dosing tubing and antidiffusion tip.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01508.220

Qt.	Order no.	Description
1 PCS	6.2803.010	Grease (2 g)
		Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX

