



## Aspiration tip automation

6.01543.100

Aspiration tip for automated OMNIS titrations.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery 6.01543.100

Qt.	Order no.	Description
1 PCS	6.01543.000	M8 aspiration tip

For aspirating off reagent in connection with the OMNIS Solvent Module.



1 PCS

6.01543.110

Tubing for aspiration tip automation, 3 pieces

Replacement hose for 601543100 aspiration tip.

