



Stirring propeller 13 mm ETFE

6.01900.020

Stirring propeller 13 mm ETFE with adjustable height for use with Rod Stirrer "Titrator" or Rod Stirrer "Sample Robot"

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.


These lists may be subject to change.

Scope of delivery 6.01900.020

Qt.	Order no.	Description
0 PCS	6.2803.000	Grease (1 oz.)
		Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Optional accessories

Order no.	Description	
6.01402.000	Sample beakers, clear glass, 75 mL, 30 pieces	
	Sample beakers, clear glass, 75 mL, suitable for 602041040 OMNIS sample rack, 30 pieces	
6.01403.040	Titration head 2xNS14 (P&P)	
	Titration head with 2 NS14 openings for use in Pick&Place modules of the OMNIS Sample Robot Pick&Place with 3 retracted dosing tips. Prepared for being equipped with rinsing nozzles. One aspiration tube is already retracted into the titration head.	