



OMNIS sample rack, 9 x 250 mL, (PP)

6.02041.010

OMNIS sample rack for OMNIS Sample Robot Pick&Place, suitable for 9 sample beakers. The following sample beakers can be used: 6.01400.000, 6.01400.003, 6.01400.100.

Plastic: Polypropylene (PP)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.



Scope of delivery 6.02041.010

Qt.	Order no.	Description
-----	-----------	-------------

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Optional accessories

Order no.	Description	
6.01400.000	Sample beakers, clear glass, 250 mL, 10 pieces Sample beakers, clear glass, 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.	
6.01400.003	Sample beakers, amber glass, 250 mL, 10 pieces Sample beakers, amber glass, 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.	
6.01400.100	Sample beakers, plastic (PP), 250 mL, 10 pieces Sample beakers, plastic (polypropylene, PP), 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.	