



OMNIS sample rack 50 x 6 mL

6.02041.060

OMNIS sample rack for OMNIS Sample Robot Oven, suitable for 50 sample vials (6.2419.000)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.02041.060

| Qt. | Order no. | Description |
|-------|-------------|----------------------------|
| 1 PCS | 6.02042.000 | Holder for rack label (3x) |

The holder for rack labels permits durable labeling of the rack with an unambiguous name.

