



# OMNIS sample rack 77 x 8 mm

6.02041.080

OMNIS sample rack for OMNIS Sample Robot – NIR suitable for 77 disposable vials with an optical path length of 8 mm for measurements in transmission (6.7402.240).

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 6.02041.080

Qt.	Order no.	Description
1 PCS	6.02042.000	Holder for rack label (3x)

The holder for rack labels permits durable labeling of the rack with an unambiguous name.

