



Double Pt-wire electrode for sample changer (fixed cable)

6.0340.000

Indicator electrode (diameter = 5.3 mm, fixed cable 2 m) used for automated Karl Fischer titration.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.02109.010	OMNIS adapter cable (Pol.) for 6.0340.000
	Adapter cable socket F on plug P for the 6.0340.000 polarizable KF electrode, connectable to the analog measuring module from OMNIS.