



Double Pt-wire electrode for coulometry

6.0341.100

Indicator electrode used for coulometric Karl Fischer titration.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.


Scope of delivery 6.0341.100

| Qt. | Order no. | Description |
|-------|------------|----------------|
| 1 PCS | 6.2008.040 | Storage vessel |

Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Optional accessories

| Order no. | Description | |
|------------|---|---|
| 6.2104.020 | Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F). |  |
