



OMNIS pipetting equipment

6.05002.010

Complete accessory set for converting the OMNIS Sample Robot Pick&Place into a version with pipetting options. The set can be mounted on all versions of the OMNIS Sample Robot (S,M and L).

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05002.010

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1562.090

Pipetting tubing / 10 mL

Pipetting tubing for 6.1562.120 and 6.05002.010.

