



## NH<sub>3</sub>-selective gas membrane electrode (low conc.)

6.0506.100

Combined ammonia-selective gas membrane electrode for low concentrations  $<10^{-2}$  mol/L, e.g.:

- in drinking water
- in boiler feed water
- in mineral water

If necessary, the complete membrane module can be easily exchanged. They are individually tested and certified.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

Order no.	Description
-----------	-------------

6.1255.000

Replacement Kit for NH<sub>3</sub>-selective Electrode 6.0506.100

Consisting of: replacement module (3 pieces), measuring electrolyte (50 mL)



---

6.2104.020

Electrode cable / 1 m / F

For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



---

6.2316.030

Inner filling solution to 6.0506.100 electrode

Measuring electrolyte for NH<sub>3</sub>-selective electrode (50 mL)

