



## Consumable Kit seals for the 6.02710.010 KF Dis-Cover lid

6.05700.170

Replacement seals for 6.02710.010 lid, 16 pieces

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Optional accessories

Order no.	Description
6.02710.010	KF Dis-Cover lid for OMNIS 120 mL sample beaker, 16 pieces
	Lid with seals for volumetric Karl Fischer Titration in 6.01400.200, 6.01400.300 and 6.01400.303 sample beakers in the OMNIS Sample Robot Pick&Place system.

