



## Platinum electrode tip

6.09395.024

Electrode tip made of platinum (Pt) for Autolab RDE, 3 mm electrode disc diameter, M4 thread connection. Temperature range: 0 - 40°C.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery 6.09395.024

Qt.	Order no.	Description
1 PCS	6.2802.000	Polishing set for solid-state electrodes
		Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)

