



Pt1000 temperature sensor for CVS (installation length 9 cm)

6.1110.120

Pt1000 temperature sensor (class B) for Professional CVS Instruments.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1110.120

Qt.	Order no.	Description
1 PCS	6.1236.020	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm and O-ring.



Optional accessories

Order no.	Description
6.2104.140	<p>Electrode cable / 1 m / 2 x 2 mm</p> <p>For connecting temperature sensors with Metrohm plug-in head G with 780, 781, 826, 827 pH/Ion meters and all Titrandos</p>

