



## Silver Electrode Tip

6.1204.130

Metal electrode tip, silver (Ag). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

Order no.	Description
6.2103.110	Contact pin M3, 2 mm  For 6.1241.060 Electrode shaft and 6.1204.1XX Electrode tips. With M3 thread. Used with VA equipment sets and for ELCD applications.



6.2802.000

Polishing set for solid-state electrodes

Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3  $\mu\text{m}$ )

