



## Platinum pin (Pt)

6.1247.010

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.


Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

Order no.	Description
6.1241.020	Electrode holder  Electrode holder for VA auxiliary electrodes with lock bolt. For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).



6.1241.120

Electrode holder

Electrode holder for VA auxiliary electrodes with plug contact 5.5. mm.  
For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).

