




MME measuring head for Professional VA instruments


6.1256.020

Measuring head for operation with the Multi-Mode electrode pro or the scTRACE Gold.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1256.020

Qt.	Order no.	Description
1 PCS	6.1829.030	PTFE tubing / M6 / 21 cm, tapered




1 PCS

6.1911.000

Flexible drive shaft for stirrer

10 pcs, for Professional VA instruments, including wrench.



1 PCS

6.2069.000

Park station for the measuring head of Professional VA/CVS instruments



1 PCS

6.2709.110

Stopper for the electrode opening on Professional VA/CVS instruments

Set with 1 stopper for electrode opening (e.g. for pH electrode or Pt1000 temperature sensor), 6 stoppers for openings (auxiliary solutions and/or automatic sample addition).



Optional accessories

Order no.	Description	
6.1454.130	Form seal made of FPM for 61414xxx for KF and VA instruments	
