



## Au working electrode, 3 mm

6.1257.210

Working electrode made of gold for amperometric detection, 3 mm diameter.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

Order no.	Description
6.1257.270	<p>Cu working electrode, 2 mm</p> <p>Working electrode made of copper for amperometric detection, 2 mm diameter. Alternative working electrode for the gold electrode in carbohydrate analysis.</p>

