



Heating tubing stopper

6.1446.210

Heating tube stopper for reagent replacement with an OMNIS Dosing Module, for 1 titration cell (coulometry)

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1446.210

Qt.	Order no.	Description
1 PCS	6.2713.000	PTFE sleeve SGJ14
		For fat-free ground-joint connection with the standard ground-joint SGJ 14.

