



# PEEK SGJ 29 stopper for KFC titration cell

6.1446.240

Stopper for volumetric titration in a titration cell (coulometry)


Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery 6.1446.240

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS	6.01406.010	Adsorber tube with lid
-------	-------------	------------------------

Adsorber tube with lid for OMNIS Karl Fischer titrator, volumetric



1 PCS

6.02709.010

Stopper

Plastic stoppers, suitable for 2.1001.4220 OMNIS Karl Fischer Titrator.



---

1 PCS

6.02709.030

Screw nipple

Plastic screw nipple, suitable for 6014050X0 titration vessel lids.



---

1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



---

1 PCS

6.2713.010

PTFE sleeve SGJ29

Sleeve for fat-free ground-joint connection with the standard ground-joint SGJ 29.

