



# Aluminum septum caps / 1000 pieces

6.1448.050

Including Al sealing.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Optional accessories

Order no.	Description
6.2419.000	Sample glass, 6 mL, 1000 pieces
	Sample glass for up to 6 mL sample, suitable for 6.1448.050 septum, 1,000 pieces



6.2621.110

Crimping tongs

Crimping tongs for sealing 6 mL sample vials (6.2419.000, 6.2419.007).

