



Screw cap with septum, 100 pcs.

6.1448.067

Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1448.077	Septum, 100 pcs.
	Replacement septum for 6.1448.067 screw cap, 100 pieces



6.2420.007

Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces


