



807 Dosing Unit 2 mL without accessories

6.1580.120

807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.1580.120

| Qt. | Order no. | Description |
|-------|------------|---|
| 1 PCS | 6.2803.010 | Grease (2 g) |
| | | Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX |

