



Drying flask cap

6.1602.145

Screw cap with GL 45 thread for the 6.1608.050 drying flask.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.1608.050	Drying bottle / 100 mL / GL 45



6.1608.090

Bottle / 250 mL / GL 45



6.1821.040

Filter tube

Filter tube for 6.1608.050 Drying bottle. For Rancimats and Karl Fischer
ovens.

