



Electrode holder with stand base

6.2001.140

Electrode holder for 949 pH Meter

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

Discover the 949 pH Tabletop Meter – Your compact pH measuring instrument for quick and reliable laboratory measurements.

- Complete pH measuring instrument with color display
- Precise pH, mV and ORP measurements with integrated temperature display
- Simple calibration with Metrohm, USA, and NIST buffers up to 3 points or by 2-point calibration with in-house buffers
- Clearly structured display of the calibrated values using symbols for simple read-off
- 3 stability criteria for measured value recognition for reliable results
- GLP conformity of the measured values with embedded date and time
- Overview of calibration data with monitoring function
- RS232 port enables connection of a printer.



Included in the package are the 949 pH Meter, the electrode holder, and the power supply unit for immediate use in the laboratory.

The 949 pH Meter is the perfect tabletop pH meter for quick and precise pH measurements. Ideal for laboratories, research institutes and industrial environments where pH analyses are performed regularly.

- Complete tabletop pH measuring instrument with color display
- Measures pH, mV, and ORP with temperature display for comprehensive analyses
- Simple pH calibration with Metrohm, USA, and NIST buffers up to 3 points or by user-defined 2-point calibration
- Display of the calibrated values with symbols for clear readability
- Signals measurement stability and offers the option of selecting between 3 stability criteria for reliable results
- In conformance with GLP with embedded date and time
- Display of calibration data and setting of the calibration deadline
- Solitrode with stable plastic shaft for long-lasting use



Included in the package are the 949 pH Meter, the electrode holder, the power supply unit, and the Solitrode with Pt1000 with fixed cable for simple and immediate use.