



# Sample rack 48 x 11 mm to IC Sample Center

6.2041.200

Sample rack for the IC Sample Center with 48 sample positions. For sample vials 6.2743.1xx.


Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Optional accessories

Order no.	Description
6.2743.100	PP Vials 0.7 mL to IC Sample Center, 1000 pieces 1000 sample vials.



6.2743.107 PP Vials 0.7 mL to IC Sample Center, 100 pieces  
100 sample vials.



---

6.2743.110 PP Vials 0.3 mL to IC Sample Center, 1000 pieces  
1000 sample vials.




---

6.2743.117 PP Vials 0.3 mL to IC Sample Center, 100 pieces  
100 sample vials.



---

6.2743.120 Septum to 6.2743.1XX Si/PTFE  
1000 pieces



6.2743.127 Septum to 6.2743.1XX Si/PTFE

100 pieces



---

6.7303.100 Glass vial, 1.5 mL, 100 pieces

For use with the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ) and 889 Sample Center.



---

6.7303.110 Septum for 67303100, 100 pieces

For use with the 67303100 Glass vial on the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ), and 889 Sample Center.

