



## Sample rack, 12 x 250 mL, (PVC)

6.2041.310

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 12 sample beakers. The following sample beakers can be used: 6.1432.320, 6.1432.323, 6.1453.220, 6.1453.250.

Plastic: Polyvinyl chloride (PVC)

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

Order no.	Description
-----------	-------------

6.1432.320 Sample beaker, clear glass, 250 mL, 1 piece

Sample beaker, clear glass, 250 mL, suitable for sample racks:  
6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,  
6.2068.020 and 6.2068.030, 1 piece



---

6.1432.323 Sample beaker, amber glass, 250 mL, 1 piece

Sample beaker, amber glass, 250 mL, suitable for sample racks:  
6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,  
6.2068.020, and 6.2068.030, 1 piece



---

6.1453.220 Sample beaker, plastic (PP), 200 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for sample  
racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,  
6.2068.020, and 6.2068.030, 1 piece



---

6.1453.250 Sample beaker, plastic (PP), 250 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample  
racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,  
6.2068.020, and 6.2068.030, 1 piece



6.1458.010

Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings.

With 6xSGJ 14 and 3x SGJ 9 openings.

