



# Power supply unit for 949 pH Meter

6.2166.200

Tabletop power supply unit, 100–240 V, 50–60 Hz, 5 V DC, for 949 pH Meter

Efficiency Level VI

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description
-----------	-------------

Discover the 949 pH Tabletop Meter – Your compact pH measuring instrument for quick and reliable laboratory measurements.

- Complete pH measuring instrument with color display
- Precise pH, mV and ORP measurements with integrated temperature display
- Simple calibration with Metrohm, USA, and NIST buffers up to 3 points or by 2-point calibration with in-house buffers
- Clearly structured display of the calibrated values using symbols for simple read-off
- 3 stability criteria for measured value recognition for reliable results
- GLP conformity of the measured values with embedded date and time
- Overview of calibration data with monitoring function
- RS232 port enables connection of a printer.



Included in the package are the 949 pH Meter, the electrode holder, and the power supply unit for immediate use in the laboratory.

---