



# Buffer solution pH 9.00 (30 x 30 mL)

6.2307.220

Buffer solution pH 9.00 (25 °C) in single use sachets, colorless

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Optional accessories

Order no.	Description
6.2307.100	Buffer solution pH 4 (500 mL)
	Ready-to-use buffer solution pH = 4 (25 °C), colored, with closure seal



6.2307.110 Buffer solution pH 7 (500 mL)  
Ready-to-use buffer solution pH = 7 (25 °C), colored, with closure seal



---

6.2307.120 Buffer solution pH 9 (500 mL)  
Ready-to-use buffer solution pH = 9 (25 °C), colorless, with closure seal



---

6.2307.200 Buffer solution pH 4.00 (30 x 30 mL)  
Buffer solution pH 4.00 (25 °C) in single use sachets, colorless



---

6.2307.210 Buffer solution pH 7.00 (30 x 30 mL)  
Buffer solution pH 7.00 (25 °C) in single use sachets, colorless



6.2307.230

Buffer solutions pH 4, 7 and 9

Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL



---

6.2325.000

pHit kit

Maintenance kit for pH electrodes

The kit contains:

- 50 mL cleaning solution
- 50 mL 3M KCl solution
- 50 mL storage solution
- 2 Storage vessels
- Instructions for use

