



## Pick-up needle MF

6.2624.100

Coated steel needle for metalfree pick-up sampling as well as for Eppendorf sample vials.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

Order no.	Description
6.2629.000	Screw nut to needle holder IC On Sample Processors



6.2833.010

Needle/tube holder 1/16 in.

For IC Sample Processors. For 1/16 in. tubes

